

Comparison sheet - Photizo Home Care units



Product Name	Vetcare	Pain Relief	Sport	Skin Care	Blush
Diode Technology	Light Emitting Diode (LED)	Light Emitting Diode (LED)	Light Emitting Diode (LED)	Light Emitting Diode (LED)	Light Emitting Diode (LED)
CE Certification	CE: For Animal use	Medical CE: Chronic Pain	Medical CE: Non-Chronic Pain	Medical CE: Wound Healing	Medical CE: Beauty Application
Dose Duration	31 seconds	41 seconds	31 seconds	16 seconds	16 seconds
Power Density (Peak / Average)	340m W/cm2 194mW/cm2	340m W/cm2 194mW/cm2	340m W/cm2 194mW/cm2	340m W/cm2 194mW/cm2	340m W/cm2 194mW/cm2
Wavelengths	633nm (visible red) 850nm (invisible NIR) both delivered simultaneously	633nm (visible red) 850nm (invisible NIR) both delivered simultaneously	633nm (visible red) 850nm (invisible NIR) both delivered simultaneously	633nm (visible red) 850nm (invisible NIR) both delivered simultaneously	633nm (visible red) 850nm (invisible NIR) both delivered simultaneously
Pulse / Continuous	30% pulsed - 70% continuous	30% pulsed - 70% continuous	30% pulsed - 70% continuous	30% pulsed - 70% continuous	30% pulsed - 70% continuous
Frequencies	5 unique frequencies	10 unique frequencies	5 unique frequencies	7 unique frequencies	7 unique frequencies
Dose delivered	6 joules/cm2	8 joules/cm2	6 joules/cm2	3 joules/cm2	3 joules/cm2
Perform up to 75.000 hours	✓	✓	✓	✓	✓
Charging Time	3-4 hrs	3-4 hrs	3-4 hrs	3-4 hrs	3-4 hrs
Mains Power	100 – 240VAC	100 – 240VAC	100 – 240VAC	100 – 240VAC	100 – 240VAC
Weight	280g	280g	280g	280g	280g
Where to buy	Click here	Click here	Click here	Click here	Click here

For more information on how to become a reseller, or how you should treat your pain, injury or wounds please don't hesitate to contact us on: post@photizo.info